

# Giraffe Shuttle

The power to provide continuity of care

The Giraffe Shuttle is a transportable power source that is an accessory to Giraffe and Panda<sup>1</sup> beds. It connects to the bed and provides electrical power to the bed and other auxiliary equipment, as well as provisions for medical gas management, required for patient care during intra-hospital transport.<sup>2</sup> The Giraffe Shuttle facilitates the mobility of these beds within a hospital building.

This reduces the potential for clinical problems associated with intra-hospital transport that result from interrupted patient thermal regulation. It also eliminates the need to transfer the infant to and from a transport incubator, reducing the potential for clinical problems associated with patient touch, handling and movement within a hospital building.

<sup>1</sup> The Giraffe Shuttle is not intended for use on Panda IWS platform.

<sup>2</sup> The Giraffe Shuttle is not intended for use outside of the hospital building or in any area where fluctuation in temperature, exposure to water, or other hazards may occur.



# Giraffe Shuttle

## The next generation of intra-hospital transport

The Giraffe Shuttle offers up to 45 minutes of electrical power to Giraffe or Panda beds, allowing you to provide continuity of care for sick babies while helping to set a new paradigm for simplified efficient infant transport within a hospital building.

- Eliminate the need to transfer a baby from bed-to-bed.**  
 The Giraffe Shuttle, in combination with the Giraffe or Panda bed, helps reduce the potential for clinical problems that can result from interrupted patient thermal regulation, patient touch, handling and movement, potentially challenging physiological stability when moving babies to and from a transport incubator.
- Customize your intra-hospital transport system.**  
 The Giraffe Shuttle provides the capability to mount various accessories and auxiliary equipment.
- Enjoy maneuverability.** The Giraffe Shuttle connects smoothly with any Giraffe or Panda bed. Designed to make it easier for you to steer the bed in tight spaces, through doorways, down hallways, and around corners.
- Achieve workflow efficiency.** In addition to all of its clinical advantages, the Giraffe Shuttle can help save time spent on transferring a baby from one bed to another and cleaning those beds.



## No more bed-to-bed transfers



The Giraffe beds combined with Giraffe Shuttle eliminate the stress of moving sick babies from one bed to another.

## One Baby. One Bed.

Consider this scenario: the next time you get notified about a high-risk delivery, just connect the Giraffe Shuttle to the Giraffe OmniBed, bring it to Labor and Delivery and admit the baby into a Giraffe OmniBed for a simplified journey from the Labor and Delivery to the NICU.

One Baby. One Bed. It's that easy. Use the same streamlined approach to move the baby from the NICU to other care areas within the same hospital building, when needed.



Giraffe Shuttle with Giraffe OmniBed

Need to leave the NICU for a scheduled procedure or diagnostic study? There's no need to disturb the baby.

Just connect the Giraffe Shuttle to the Giraffe Incubator to move the baby to the desired care area within the same hospital building.



Giraffe Shuttle with Giraffe Incubator

The Giraffe Shuttle provides power to the Giraffe Warmer, enabling clinicians to transport the baby already in the warmer from one care area to another.



Giraffe Shuttle attaching to Panda Warmer